



IDAHO DEPARTMENT OF
HEALTH & WELFARE

Test Title	<i>Campylobacter</i> spp. Identification/Confirmation of Referred Isolate Methodology: Biochemical Testing, Sequencing of 16S Ribosomal DNA
Specimen Requirements	1. Specimen type(s): actively growing pure culture on suitable medium, preferably in microaerophilic atmosphere
Sampling Materials	1. Sample Container: sealed culture tubes or plates
Procedural Notes	1. Bacteriology Test Request Form 2. CPT Code: 87077 (Culture), 87158 (Sequence Analysis)
Shipping Instructions	1. Temperature/preservative instructions: ship at room temperature in microaerophilic transport environment. 2. Package according to Biological Substance, Category B, shipping guidelines. 3. Ship to: Idaho Bureau of Laboratories ATTENTION: Bacteriology Laboratory 2220 Old Penitentiary Rd Boise, ID 83712
Reporting and Turnaround Time (TAT)	1. TAT: identification from pure culture available in 4-6 working days 2. This disease must be reported to your local health district or to the state Bureau of Communicable Disease Prevention according to the rules and regulations governing Idaho reportable diseases (IDAPA 16.02.10).